

# TASKI® Sani 4 in 1 Plus Conc.

## Washroom cleaner, descaler, disinfectant and deodoriser in SmartDose® BPR approved

### Description

Washroom cleaner, descaler, disinfectant and deodoriser for all acid resistant hard surfaces in washroom areas.

### Key properties

- Acid based formulation
- Cleans, descales, disinfects and deodorises in one step
- Effective against bacteria, enveloped viruses and yeast
- Unique, patented Odour Neutralizing Technology (O.N.T.)

### Benefits

- Excellent cleaner, works in a large spectrum of water hardness
- Reduces the number of products in the washroom area
- Excellent hygiene, reducing the risk of cross contamination
- Leaves a pleasant, fresh fragrance

### Use instructions

TASKI Sani 4 in 1 Plus is supplied in SmartDose dosing system for convenient and controlled dosing of the product into a spray bottle or bucket.

### Dosage\*

Spray application: Fill the 500 ml trigger bottle with water. Turn yellow head to trigger bottle position, pull head up, then push down fully to dispense 1x10ml (2%) into the trigger bottle. Attach the trigger and shake bottle gently.

Bucket application: Fill bucket with water. Turn yellow head to bucket position: pull head up, then push down fully. Dose 1x70ml per 3.5L water (2%).



*\*This dosage is according to optimal conditions, recommendations may vary, please consult with your Diversey representative for directions.*

### Application

Outside of healthcare environments (e.g. public areas, hotels and other institutional areas such as schools/offices etc.) the product use is approved in a 1-step cleaning process (dirty conditions).

In healthcare environments, defined where disinfection is medically indicated (e.g. doctor offices, ambulances, hospital wards and LTC), the use of the product requires a 2-step cleaning process. Remove gross soiling and pre-clean the surface before applying the product.

- Apply the diluted product by fully wetting the surface.
- Use a cloth, brush, sponge or mop for even distribution.
- Leave the surface wet for 5 minutes against bacteria and enveloped viruses and 15 minutes against yeast.
- Wipe dry or let the surface air dry.
- Rinsing of the surface is not required.



# TASKI® Sani 4 in 1 Plus Conc.

## Technical data

Appearance: Clear red liquid  
pH value (neat): < 2  
pH value (in use): < 2  
Relative density (20°C): ≈ 1.09

*The above data is typical of normal production and should not be taken as a specification.*

## Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diversey.com. Store in original closed containers away from extremes of temperature. Only for professional users / specialists. Use disinfectants safely. Always read the label and product information before use and follow local regulations and advisory.

## Product compatibility

Do not mix with other products. Do not use on water or acid sensitive surfaces. Sufficiently pre-wet grout with water before use. Test material compatibility in a small, inconspicuous place before use.

## Microbiological data

EN Test Method	Soiling Conditions	Temp. °C	Contact Time	Test organisms	Pass Concentration % v/v
EN 1276	3.0 g/L BSA	20	5 min	Staphylococcus aureus Enterococcus hirae Escherichia coli Pseudomonas aeruginosa	1.5%
EN 1650	3.0 g/L BSA	20	15 min	Candida albicans	2%
EN 13697	3.0 g/L BSA	24	5 min	Staphylococcus aureus Enterococcus hirae Escherichia coli Pseudomonas aeruginosa Candida albicans	2%
EN 14476	3.0 g/L BSA + blood	20	5 min	Modified Vaccinia virus Ankara	1%
EN 16777	3.0 g/L BSA + blood	25	5 min	Modified Vaccinia virus Ankara	1%
EN 13727	0.3 g/L BSA	20	1 min	Staphylococcus aureus Enterococcus hirae Pseudomonas aeruginosa	1%
EN 13624	0.3 g/L BSA	20	5 min	Candida albicans	1%
EN 13697	0.3 g/L BSA	20	5 min	Staphylococcus aureus Enterococcus hirae Escherichia coli Pseudomonas aeruginosa	1%
EN13697	0.3 g/L BSA	20	5 min	Candida albicans	1%

## Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EU Detergent Regulation, EC 648/2004.

## Available pack sizes

TASKI Sani 4 in 1 Plus Conc is available in 1.4L SmartDose. Empty spray bottles and corresponding sprayers are also available.